



DISPOSABLE ECG ELECTRODE



SF01



SF05



SF08 / BA08



SF11



SN08



ST31



SC10

General Monitoring

Foam Stainless Steel Snap

Medical grade material is impervious to fluids, high tack medical grade adhesive is ideal for long term monitoring

Recommended usage

ICU/CCU, Surgery, Emergency

Item#	Description	Packaging	Case Dimension
SF11	Child, Solid Gel, Ø30mm	50/pk, 40pk/cs	48×36×28cm
SF03	Adult, Solid Gel, 50×55mm	50/pk, 40pk/cs	48×36×28cm
SF05	Adult, Solid Gel, 42×50mm	50/pk, 40pk/cs	48×36×28cm
SF06	Adult, Solid Gel, Ø55mm	50/pk, 40pk/cs	48×36×28cm
SF08	Adult, Solid Gel, Ø50mm	50/pk, 40pk/cs	48×36×28cm
BA08	With 4mm banana plug, Ø50mm	50/pk, 40pk/cs	43×55×28cm
SF07	Neonatal, Solid Gel, Ø25.4mm	50/pk, 40pk/cs	40×36×28cm

Cloth Stainless Steel Snap

Soft and breathable cloth material, particularly designed for sensitive skin.

Recommended usage

ICU/CCU, Emergency, Long term

Item#	Description	Packaging	Case Dimension
SN08	Cloth, Solid Gel, Ø50mm	50/pk, 40pk/cs	40×36×28cm

Clear Tape Stainless Steel Snap

Breathable Tape allows caregiver to monitor the patient's skin directly under the electrode without having to remove it. Low profile design promotes patient comfort.

Recommended usage

ICU/CCU, Surgery, Emergency, Long term, Holter

Item#	Description	Packaging	Case Dimension
ST31	Breathable clear tape, Solid Gel, Ø47X43mm	50/pk, 40pk/cs	40×36×28cm

Radiotranslucent

Foam Radiolucent Snap

Recommended usage

Surgery, Emergency, Cath Lab, Radiology, X-ray, CT and MRI applications.

Item#	Description	Packaging	Case Dimension
SC10	Adult, Solid Gel, Ø42X37mm	50/pk, 40pk/cs	48×36×28cm

Short Monitoring, Long Monitoring, Stress Test, Holter, emergency procedures, Radio Translucent.

ALL THE ECG ELECTRODES ARE NON-IRRITATING; NON-SENSITIZING; NON-CYTOTOXIC; LATEX-FREE AND PVC-free.



Stress/Diaphoretic

Intended for use in applications where high tack adhesive is needed.

Recommended usage

ICU/CCU, Surgery, Stress, Holter, Emergency

Item#	Description	Packaging	Case Dimension
SF31	Adult, Solid Gel, 52x42mm	50/pk, 40pk/cs	48x36x28cm



SF31

Wet Gel

Conductive wet gel allows for a quick, accurate EKG signal with high tack adhesive.

Recommended usage

ICU/CCU, Surgery, Stress, Holter, EMS, Emergency

Item#	Description	Packaging	Case Dimension
WF32	Adult, Wet Gel, 52x47mm	50/pk, 40pk/cs	48x36x28cm



WF32

Prewired Neonatal

Manufactured with a conductive adhesive hydrogel and breathable, stretchable, gentle materials, ideal for neonatal patients.

Recommended usage

NICU, PICU

Item#	Description	Packaging	Case Dimension
PW01	Neonatal, 22x22mm	3/pk, 10pk/bx, 12bx/cs	48x22x22cm



PW01

Repositionable

Manufactured with a conductive hydrogel, allowing the caregiver to reposition at least one time.

Recommended usage

ICU/CCU, Surgery, Emergency, Cath Lab, Radiology

Item#	Description	Packaging	Case Dimension
RS01	Adult, Solid Gel, 32x40mm	3/pk, 10pk/bx, 12bx/cs	48x22x22cm



RS01

Diagnostic (Resting) Disposable ECG Electrodes.

With conductive and adhesive hydrogel, reduce potential for cross-contamination, particularly intended for short term monitoring (Emergency Room)

Recommended usage

Cardiology, resting ECG

Item#	Description	Packaging	Case Dimension
TAB03	20x34mm	10/pk, 50pk/bx, 10bx/cs	30x19x23cm
TAB04	23x34mm	10/pk, 50pk/bx, 10bx/cs	30x19x23cm



TAB03



TAB04

Short Monitoring, Long Monitoring, Stress Test, Holter, emergency procedures, Radio Translucent.

ALL THE ECG ELECTRODES ARE NON-IRRITATING; NON-SENSITIZING; NON-CYTOTOXIC; LATEX-FREE AND PVC-free.

Imported by :  **ASIAMEDIC**
Hospital Supply Sdn Bhd
1154902-T

T : +6012 - 5050 158 / +6017 - 4222 236
F : +604 - 6192 699

E : info@asiamedic.com.my
W : www.asiamedic.com.my