

Dual Lumen Catheters and kits



Offering optimal balance of flexibility, anti-kinking, durability, size and flow.

Radiopaque biocompatible tecoflex polyurethane sterile catheter, different lengths and sized that determine the suitability of blood flow and patient age.

- Heavy duty length marked guide wire with advancer.
- Long tapered tip dilators with stopper to avoid over insertion.
- Scalpel # 11 for small snipping.
- Guiding Y Connector “GY” (option).
- Nitinol Guide wire.

Ordering Information:

Fr Size	Length	Subclavian	Jugular	Pre Curved
		KIT CODE	KIT CODE	KIT CODE
6	8 Cm	DLC-0608-K		
7	10 Cm	DLC-0710-K	DLC-0710-JK	
7	13 Cm	DLC-0713-K	DLC-0713-JK	
8	10 Cm	DLC-0810-K	DLC-0810-JK	
8	13 Cm	DLC-0813-K	DLC-0813-JK	
9	10 Cm	DLC-0910-K	DLC-0910-JK	
9	13 Cm	DLC-0913-K	DLC-0913-JK	DLC-0913-PCK
9	15 Cm	DLC-0915-K	DLC-0915-JK	DLC-0915-PCK
12	12 Cm	DLC-1212-K	DLC-1212-JK	DLC-1212-PCK
12	15 Cm	DLC-1215-K	DLC-1215-JK	DLC-1215-PCK
12	17 Cm	DLC-1217-K	DLC-1217-JK	DLC-1217-PCK
12	20 Cm	DLC-1220-K	DLC-1220-JK	DLC-1220-PCK
14	15 Cm	DLC-1415-K	DLC-1415-JK	DLC-1415-PCK
14	17 Cm	DLC-1417-K	DLC-1417-JK	DLC-1417-PCK
14	20 Cm	DLC-1420-K	DLC-1420-JK	DLC-1420-PCK