

## Triple Lumen Catheters kits



Offering optimal balance of flexibility, anti-kinking, durability, size and flow.

Radiopaque biocompatible tecoflex polyurethane sterile catheter, different lengths and sized that determine the suitability of blood flow and patient age.

1. Heavy duty length marked guide wire with advancer.
2. Long tapered tip dilators with stopper to avoid over insertion.
3. Scalpel # 11 for small snipping.
4. Guiding Y Connector “GY” (option).
5. Nitinol Guide wire “N” (option).

### Ordering Information:

Fr Size	Length	Subclavian	Jugular	Pre Curved
		KIT CODE	KIT CODE	KIT CODE
13	15 cm	TLC-1315-K	TLC-1315-KJ	TLC-1315-PCK
13	17 cm	TLC-1317-K	TLC-1317-KJ	TLC-1317-PCK
13	20 cm	TLC-1320-K	TLC-1320-KJ	TLC-1320-PCK
14	15 cm	TLC-1415-K	TLC-1415-KJ	TLC-1415-PCK
14	17 cm	TLC-1417-K	TLC-1417-KJ	TLC-1417-PCK
14	20 cm	TLC-1420-K	TLC-1420-KJ	TLC-1420-PCK