

Introduction

- The AutoFuser® pump is intended for continuous and/or intermittent infusion of medication for general infusion use, including antibiotic, chemotherapy and pain management therapies
- Included strap clip and carrying pouch allow for ambulatory flexibility, regardless of catheter site location
- Unlike some elastomeric pumps, the AutoFuser® is easy to prepare as the pump can lie flat during the filling process for simple syringe or filling device transfer
- The recessed fill port minimizes sterility issues (Within minutes, the pump can be filled, primed, and ready for patient use)



Advantages

- Portable and disposable silicone infusion pump
- Pain relief without the side effects of narcotics
- Optional PCA bolus volumes (0.5, 1.0, 2.0 and 5.0ml to address breakthrough pain)
- Unique shape provides easy filling and patient comfort
- Bolus button is patient controlled for breakthrough pain control
- Reduce hospital stays and improve recovery time
- Latex free and DEHP free

Specification

AutoFuser	
Type	Continuous type, bolus type, continuous & bolus type
Volume size	60ml, 100ml, 150ml, 275ml, 400ml and 550ml reservoirs available
Basal rates	0.5 - 10.0ml/h (0.5, 1.0, 2.0, 4.0, 5.0, 8.0, 10.0ml/h)
Bolus volume	0.5, 1.0, 2.0 and 5.0ml to address breakthrough pain
Lockout times	10, 15, 30 and 60min.