

## Selective Flow Rate Infuser

### Introduction

- The AutoSelector® is a disposable balloon pump with an integrated flow regulator that allows for flexible programming
- Multiple choices in flow rates and PCA options allow you to customize your patient's pain protocol
- Reservoir design significantly reduces pressure needed to fill balloon
- Local anesthesia pain management promotes (Increased patient mobility, Faster recovery, Fewer readmissions, Decreased narcotic usage, Lower infection rates, Greater patient satisfaction)



### Advantages

- Portable and disposable silicone pumps
- Two types of pumps (A Type 0.5ml/h : 7.0ml/h, B Type 1.0ml/h : 14.0ml/h)
- The AutoSelector® light weight and small size make them more comfortable for the patient to carry when ambulatory
- Safety lock with removable key and lockable cover prevents patient tampering
- 8 incremental flow rates
- Comfortable to use, ensuring patient compliance
- Front fill port accommodates easy filling
- Available in a various volume size (60ml, 100ml, 150ml, 275ml, 400ml, 550ml)
- Rounded, low-profile housing provides
- Latex free and DEHP free

### Specification

AutoSelector	
Type	Continuous model, Continuous & bolus Model
Volume size	60ml, 100ml, 150ml, 275ml, 400ml, 550ml
Flow rate	A type : 0.5ml/h ~ 7.0ml/h, B type : 1.0ml/h ~ 14.0ml/h
Bolus volume	0.5, 1.0, 2.0, 5.0ml
Lockout times	10, 15, 30 and 60min.