

## Introduction

- The Autofuser S® is ideal for the flexible delivery requirements of Postoperative pain management, Cancer pain management, Chemotherapy, and Antibiotic therapy
- Portable, compact, lightweight and discreet fits into a pocket or bag maximizing patient compliance and convenience



## Advantages

- The Bolus button of Autofuser S® can help you to save priming time by pressing bolus in advance
- Optimized fill port to smoothly fill medication while minimizing resistance
- Latex and DEHP free for clinician and patient safety
- Protective housing for a crush-resistant, non-disruptive medication infusion
- Medical grade non-kinking PVC
- Clamp and carrying pouch for convenient securement
- Air filter for removing the residual air in and tubing

## Specification

AutoFuser S	
Flow rate	2.0 - 4.0ml/h
Bolus dose	0.5, 1, 2, 5ml
Lockout time	1, 15, 30, 60min.
Size	100ml, 150ml, 275ml, 400ml, 550ml