

# GT400 Disposable Automatic Biopsy Instrument

GT400 can be easily operated with one hand and is an ideal instrument for ultrasound guided biopsy procedures.

## Features

- A 22mm throw with 18mm sample notch harvests quality soft tissue specimens.
- Two firing triggers, rear and side trigger accommodate different preference.
- Built-in Safety prevents accidental firing to improve safety.
- Color coded body for easy and accurate needle gauge identification.
- Ergonomic design and lightweight
- Compatible with coaxial needles.



	Gauge	Length 10cm	Length 15cm	Length 20cm	Length 25cm
GT400 Disposable Automatic Biopsy Instrument	14	G-T211410	G-T211415	G-T211420	-
	16	H-171610	H-171615	H-171620	-
	18	G-T211810	G-T211815	G-T211820	B-1211800
	20	D-T212010	D-T212015	D-T212020	-

# GTC Coaxial Needle for GT400 and GT20

Specially designed for use with GMT's GT400 and GT20 biopsy Instruments.

## Features

- Engineered compatibility with GT400 and GT20 biopsy instruments facilitates efficiency and accuracy
- Lightweight and ultra thin wall construction
- Creates a clear path to work through when performing multiple biopsies in the same area
- Sized just one gauge larger than corresponding core biopsy needles
- Color-coded depth stopper conveniently matches the compatible GT400 and GT20 biopsy instruments.



	Gauge	Length 10cm	Length 15cm	Length 20cm
GTC	14	B-1121110	B-T121415	B-T121420
	16	D-T121610	B-1121615	B-1121620
	18	B-T121810	B-T121815	B-T121820
	20	B-T122010	B-1122015	B-1122020
	21	B-1122110	B-T122115	B-T122120

Imported by :



**ASIAMEDIC**  
Hospital Supply Sdn Bhd  
1154902-T

T : +6012 - 5050 158 / +6017 - 4222 236  
F : +604 - 6192 699  
E : info@asiamedic.com.my  
W : www.asiamedic.com.my